

Clear Vial™

5 mL, 10 mL, 20 mL, 50 mL, 100 mL



Ref. no.	1875	1876	1877	1879	1880
Volume (mL)	5	10	20	50	100
Material	COC	COC	COC	COC	COC
Colour	Natural	Natural	Natural	Natural	Natural
Dimension A (mm)	40	45	55	73	94.5
Dimension B (mm)	Ø22	Ø24	Ø30	Ø42.5	Ø51.6
Dimension C (mm)	Ø20	Ø20	Ø20	Ø20	Ø20
Weight (gr)	2.9	3.5	5.1	12.5	19.5
Brimful Volume (mL)	10	14	27	70	134
Packing (pcs/pack unit)	230	162	135	55	36
Packing (pcs/carton)	3,680	2,592	1,620	770	432
Volume of carton (m³)	0.094	0.094	0.094	0.144	0.144
Carton dimension LxWxH (mm)	600x400x391	600x400x391	600x400x391	800x400x450	800x400x450
Carton weight filled (kg)	11.5	9.9	9.0	10.4	9.2
Cartons/EUR pallet	12	12	12	9	9
Quantity/EUR pallet	44,160	31,104	19,440	6,930	3,888
Sterilization	Irradiation/Autoclaving				
Rubber disc or stopper - Cap	The vials fit rubber plugs according to ISO 8362-2 and aluminium caps according to ISO 8362-3/6				

2024-09

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.