


Vials for Injectables Polypropylene (PP)

3 mL, 5 mL, 10 mL, 20 mL, 25 mL



Ref. no.	1834	1835	1836	1837	1838
Volume (mL)	3	5	10	20	25
Material	PP	PP	PP	PP	PP
Colour	Natural / Amber	Natural / Amber	Natural / Amber	Natural / Amber	Natural / Amber
Dimension A (mm)	25.9	37.9	45.6	48.6	53.5
Dimension B (mm)	14	26	32	29	38
Dimension C (mm)	Ø22	Ø22	Ø26	Ø32	Ø33
Neck opening (mm)	Ø12.6	Ø12.6	Ø12.6	Ø12.6	Ø12.6
Weight (gr)	2.0	2.8	3.9	5.5	6.3
Brimful Volume (mL)	5	9	15	24	30
Packing (pcs/pack unit)	230	230	162	98	85
Packing (pcs/carton)	5,760	4,140	2,592	1,372	1,105
Volume of carton (m ³)	0.094	0.094	0.094	0.094	0.094
Carton dimension LxWxH (mm)	600x400x391	600x400x391	600x400x391	600x400x391	600x400x391
Carton weight filled (kg)	12.5	12.3	11.2	8.5	7.9
Cartons/EUR pallet	20	20	20	20	20
Quantity/EUR pallet	115,000	82,800	51,840	27,440	22,100
Sterilization	Autoclavable				
Rubber disc or stopper - Cap	The vials fit rubber plugs according to ISO 8362-2 and aluminium caps according to ISO 8362-3/6				

2024-09

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.