


# Intra-Con Single Neck IV Container

100 mL



<b>Ref. no.</b>	<b>1923</b>
<b>Volume (mL)</b>	100
<b>Material</b>	PP
<b>Additive port</b>	No
<b>Colour</b>	Natural
<b>Dimension A (mm)</b>	213.8
<b>Dimension B (mm)</b>	58
<b>Dimension C (mm)</b>	27
<b>Neck opening (mm)</b>	Ø20
<b>Weight (gr)</b>	11.2
<b>Packing (pcs/carton)</b>	552
<b>Volume of carton (m<sup>3</sup>)</b>	0.094
<b>Carton dimension LxWxH (mm)</b>	600x400x391
<b>Carton weight filled (kg)</b>	7.1
<b>Cartons/EUR pallet</b>	20
<b>Quantity/EUR pallet</b>	11,040
<b>Sterilization</b>	Autoclavable

2024-09

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.